



Home Requirement Booklet

User Information Booklet



About this booklet

Before the team can deliver a MyPAM device to your home, you will need to confirm that your space meets the following product requirements.

We ask you to consider these requirements when deciding on a suitable space for the MyPAM device. The space you choose will be the dedicated space for the device for the duration of your MyPAM rehabilitation plan.

This booklet is intended for service users who have enrolled in the MyPAM service or clinical trial. This booklet is to be given to service users before receiving the MyPAM device.

This booklet includes:

- MyPAM product footprint
- MyPAM home requirements

What is MyPAM?

MyPAM is a home-based robotic rehabilitation device, designed to deliver safe and effective upper-limb rehabilitation to those who have poor arm function or weakness as a result of stroke.

The MyPAM device consists of a powered 'arm' which uses motors to guide the user through a range of exercises.

The exercises are delivered through a series of on-screen computer games, which have been designed to increase engagement with physiotherapy.

The exercises aim to reflect those of traditional physiotherapy, helping users to regain movement and strength in their affected arm.

PRODUCT DIMENSIONS

Product Dimensions

H: 1500 (1200)mm W: 1500 (850)mm D: 1100 (700)mm

Workspace Dimensions

Height: 1500mm Width: 850mm Depth: 1300mm

Max Arm Extension


750mm

Product Weight

30kg



SIDE VIEW



PRODUCT DIMENSIONS

Product Dimensions

H: 1500 (1200)mm W: 1500 (850)mm D: 1100 (700)mm

Workspace Dimensions

Height: 1500mm Width: 850mm Depth: 1300mm

Max Arm Extension


750mm

Product Weight

30kg



TOP VIEW



PRODUCT DIMENSIONS

Product Dimensions

H: 1500 (1200)mm W: 1500 (850)mm D: 1100 (700)mm

Workspace Dimensions

Height: 1500mm Width: 850mm Depth: 1300mm

Max Arm Extension


750mm

Product Weight

30kg



AVERAGE WORK PLANE



HOME REQUIREMENTS

Before having a MyPAM device in your home, please consider the following requirements


You will need to confirm that your space meets the requirements, before the MyPAM device can be installed



1 Enough floor space


The MyPAM device is about the size of an armchair.

You will need a total floor space of **1.5 by 1.3 meters** to use the device safely.



2 Access to mains electricity

Please choose a space with easy access to a plug socket.



3 Chair with back support

You will need **your own chair** with a good back support and fixed position.

This will be the allocated chair for using MyPAM.

4 Access to internet connection

The device uses the internet to securely send progress data to your therapist.

The MyPAM team will provide a Wi-Fi hub, but you still need internet access.

HOME REQUIREMENTS

Before having a MyPAM device in your home, please consider the following requirements

You will need to confirm that your space meets the requirements, before the MyPAM device can be installed



Requirement Checklist

Don't worry if you do not have access to all of the requirements. Contact a member of the MyPAM service team to let them know. We will try and accommodate as best as possible.

What you will need

(please tick to confirm that you have the following requirements before delivery)

- Minimum floor space of **1.5 by 1.3 meters**
- A chair** with good back support *(this will be the 'designated' chair for the duration of the MyPAM service)*
- A near-by plug socket
- Internet access and ethernet port.

What we will provide

(Please tick to confirm that you have received the following items on delivery)

- A MyPAM Device with Screen
- Handle and Wrist rest
- A MyPAM Wi-Fi hub and ethernet cable
- Extension cable

Please check that you receive all the above items on delivery.



For more information and to find out how to get involved in the research, get in touch:



hello@robotic.rehab



www.robotic.rehab